



IHE Integration Statement			
Vendor	System Name	System Version	Date
Carl Zeiss Meditec	HFA II-i	5.1	2011, Oct 10
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profile Implemented	Actor Implemented	Integration Profile Option Implemented	
Basic Eye Care Workflow	Acquisition Modality	None	
Internet address for vendors IHE information: <a href="http://www.meditec.zeiss.com/dcs">http://www.meditec.zeiss.com/dcs</a>			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	<a href="http://www.meditec.zeiss.com/dcs">http://www.meditec.zeiss.com/dcs</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	