

## The cold light for the small budget

## KL 200 – The new standard for stereomicroscopes from Carl Zeiss

The wide variety of applications of stereomicroscopes makes an illumination system necessary which is matched to the applications concerned. It is important to illuminate the object not only homogeneously, but also with optimum protection, i.e. without heating it.


With the **new Schott KL 200 cold-light source**, Zeiss is rounding off its line of cold-light sources in the lower price category for stereomicroscopes.

**Its main features are:**


- **small, compact, esthetic**
- **8 V 20 W lamp, easy-to-change**
- **brightness in 3 steps**
- **excellent price-performance ratio**

Due to the versatile accessories (can be used with the KL 750), the illumination can now be adapted very effectively to the applications also in this price category. This provides the user with an additional benefit.

The KL 200 will replace the 6 V 10 W halogen illumination (485001 9804). The L stand (455108) will remain unchanged for the time being.



# System-Diagramm



For further details, please contact:



**Carl Zeiss  
Mikroskopie**  
D-07740 Jena  
Telefon: (0 36 41) 64-16 16  
Telefax: (0 36 41) 64-31 44  
Internet: mikro@zeiss.de  
<http://www.zeiss.de>